CERICHEM[®]
B I O P H A R M

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Packaging
Cod/SKU HP2100
Fiacone da 1000 ml/Bottle 1000 ml

DISINFECTANT Multi-enzymatic DETERGENT

DESCRIPTION

Ready to use plurienzymatic disinfectant for the decontamination of critical and semicritical hospital surgical instruments. To use before the transport and reprocessing (manual and machine). Also suitable for the treatment of not sumersible medical instruments.

COMPOSITION

Enzymatic mixture, Cationic surfactants, Nonionic surfactants.

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METHOD OF USE

Directly nebulize on devices' surfaces to be treated using the spray dispenser. Leave on for 15 minutes.

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Benzalkonium chloride based DISINFECTANT

DESCRIPTION

Concentrated disinfectant specific for the decontamination of critical and semicritical hospital surgical instruments, medical devices.

COMPOSITION

Cationic surfactants, nonionic surfactants, o-phenylphenol.

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METHOD OF USE

To clean the instruments dilute to 2% and immerse the instruments for a contact time of 15 minutes, then rinse.

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METHOD OF USE

Ready to use solution. Pour the content and immerse the instruments for a contact time of 15 minutes. Then rinse.

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Multi-enzymatic DETERGENT

DESCRIPTION

Concentrated multi-enzymatic detergent specific for the cleansing of critical and semicritical hospital surgical instruments, medical devices and endoscopy material.

COMPOSITION

Enzymatic mixture (5 enzymes), Cationic surfactants, Nonionic surfactants.

METHOD OF USE

For a cleaning activity, dilute to 0,5% and immerse the instruments. Rinse.

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**HIGH LEVEL DISINFECTANT
based on peracetic acid**

DESCRIPTION

PATHOSTER® 0,35% is an innovative high level disinfectant solution, stabilized, based on 0,35% peracetic acid and hydrogen peroxyde. Suitable for the disinfection of surgical instruments, endoscopy material, invasive medical devices in high biological risk environments (biomedical analysis laboratories, biological laboratories, hospitals, operating theatres, healthcare environments).

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METHOD OF USE

The device can be used:

1. by immersion;
2. nebulized;
3. in endoscope washers.

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Keywords: peracetic acid; highly pathogen agents; sterilization

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ABSTRACT

The mainstays of biological risk prevention and control of infectious disease are the disinfection and the sterilization. It is not easy to find products with ideal characteristics for different pathogen agents and situations. We evaluated the antimicrobial effectiveness and stability of two formulations of a commercial product of stabilized peracetic acid (Pathoster® 0,35% and Pathoster® 0,5%) against spores of *Bacillus anthracis* and *Bacillus cereus* and vegetative forms of *Yersinia pestis*, *Burkholderia mallei*, *Burkholderia pseudomallei*, *Francisella tularensis*, *Brucella abortus* and *Brucella melitensis*. Both products inactivated spores of *Bacillus anthracis*, *Bacillus cereus* and vegetative forms of the others bacteria after 5 minutes of exposition at room temperature for at least 2 weeks since reconstruction of the products. Neither formulation induced toxicity to humans.

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Bibliografia

1. Fretzschel H et al. 2013; Decontamination of high-risk animal and zoonotic pathogens. *Biosecur. Bioterror.* (12)pp: 102-114
2. Leggett MJ et al. 2015; Resistance to and killing by the sporicidal microbicide peracetic acid. *J Antimicrob Chemother.* 70(3)pp: 773-779
3. Tsumi M et al. 1973; Antimicrobial action of peracetic acid. *J Food Hyg Soc.* (15)pp: 116-120
4. Freer PC, Novy FG. 1902; On the formation, decomposition and germicidal action of benzoylperoxy and diacetyl peroxides. *Am Chem J.* (27)pp: 161-193
5. Block SS. 1991; Disinfection, sterilization and preservation (4 ed.). *Lee & Febiger Publ., Philadelphia.*

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**DISINFECTANT based
on peracetic acid**

DESCRIPTION

High level disinfectant (sporicidal) to activate, stabilized, based on peracetic acid 0.19% and hydrogen peroxide.

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METHOD OF USE

The device can be used (1) by immersion and (2) in endoscope washers.

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REVITALIZING for stainless steel instruments

DESCRIPTION

Revitalizing concentrated solution, descaling, efficient and fast. Removes the mineral spots that may occur after steam sterilization processes. **Lubricates the instruments.**

COMPOSITION

Phosphoric acid, nonionic surfactants.


METHOD OF USE

Mix RUGILUX to 12% in water. Put the clean instruments at the bottom of the tray. Immerse them in the solution. Remove the instruments from the solution and rinse.

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Strongly alkaline DETERGENT

DESCRIPTION

Alkaline Detergent solution for medical devices' treatment in automatic systems. Concentrated liquid.

COMPOSITION

Sequestering agents, alkaline agents. Phosphate free.

METHOD OF USE

Used in cleaning equipments 2-5 ml/l. A neutralizing phase by using H-MATIC NPC or H-MATIC NP is required.

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Slightly alkaline DETERGENT

DESCRIPTION

Concentrated detergent solution, enzymatic, slightly alkaline, specific for surgical instruments, endoscopes and medical devices.

COMPOSITION

Anionic surfactants, non ionic surfactants, enzymes, alkalinity donors.

METHOD OF USE

Used in cleaning equipments 2-5 ml/l. Does not require neutralization.

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Multi-enzymatic DETERGENT

DESCRIPTION

Concentrated plurienzymatic solution, detergent, specific for surgical instruments and medical devices. **The solution doesn't have corrosive features.**

COMPOSITION

Anionic surfactants, non ionic surfactants, enzymes.

METHOD OF USE

Used in cleaning equipments 3-10 ml/l. Does not require neutralization.

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NEUTRALIZING based on inorganic acids

DESCRIPTION

Neutralizing detergent solution based on phosphoric acid for the treatment of medical instruments in automatic systems after the cleaning phase with alkaline detergents.

COMPOSITION

Phosphoric acid, corrosion inhibitors.

METHOD OF USE

Measure between 1-3 ml/l or through a door machine incorporated dispenser.

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DETERGENT NEUTRALIZING based on organic agents

DESCRIPTION

Neutralizing detergent solution based on organic acid for the treatment of medical instruments in automatic systems after the cleaning phase with alkaline detergents.

COMPOSITION

Organic acids, sequestrants.

METHOD OF USE

Measure between 1-3 ml/l or through a door machine incorporated dispenser.

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High level DISINFECTANT

DESCRIPTION

High level disinfectant, decontaminant, based on organic acids and peroxides for automatic systems.

COMPOSITION

Peracetic acid, hydrogen peroxide.

MICROBIOLOGICAL ACTIVITY

Bactericidal: (EN 13727 – EN 14561), Fungicidal: (EN 13624 – EN 14562), Tuberculocidal: (EN 14348 – EN 14563), Virucidal: (EN 14476), Sporocidal: (EN13704).

METHOD OF USE

In order to have a high level disinfection, dose 1% of product (10ml in 1L of water) and leave on for 15 minutes.

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High level DISINFECTANT

DESCRIPTION

High level disinfectant, decontaminant, based on organic acids and peroxides for automatic systems.

COMPOSITION

Peracetic acid, hydrogen peroxide.

MICROBIOLOGICAL ACTIVITY

Bactericidal: (EN 13727 – EN 14561), Fungicidal: (EN 13624 – EN 14562), Tuberculocidal: (EN 14348 – EN 14563), Virucidal: (EN 14476), Sporocidal: (EN13704).

METHOD OF USE

Dose 3% of product (30 ml in 1L of water) for a high level disinfection. Leave on for 15 min.

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Descaler and rinse activator

DESCRIPTION

Descaler detergent and brightener, formulated with slightly foaming surfactants. It could be used as a liquid rinsing agent with lubricant effect and drying activator to use in disinfectors. Also suitable as descaler for pans washers.

COMPOSITION

Nonionic surfactants, organic acids, sequestrants.

METHOD OF USE

Measure between 0,5-1,5 ml/l or through a door machine incorporated dispenser.

H-Matic LUB

- ✓ Idoneo per strumentario chirurgico articolato
Suitable for articulated surgical instruments
- ✓ Non contiene siliconi
It doesn't contain silicones

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LUBRICANT for surgical instruments

DESCRIPTION

Lubricant for surgical instruments to use in the last rinse in automatic systems.

COMPOSITION

Aliphatic hydrocarbons, nonionic surfactants.

METHOD OF USE

Measure 1,2 ml/l of lubricant H-MATIC LUB through the appropriate dosing pump integrated in the machine.

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DISINFECTANT based on peracetic acid

DESCRIPTION

Concentrated high level disinfectant solution (sporicidal), twocomponent, based on peracetic acid and hydrogen peroxide, suitable for washing machines-disinfectors.

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
METHOD OF USE

Measure 30 ml of A solution + 30 ml of B solution for each liter of water, in order to have a peracetic acid concentration of 0,15% (1500 ppm).

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PROCESS ♦ DET + ACT + F

H-MATIC DET

is a detergent for the automatic treatment of medical thermolabile devices formulated with non-foagenic surfactants. It's used for the automatized treatment of flexible endoscopes and other medical devices present in healt facilities or biologic labs.

H-MATIC ACT

is activator for process based on peroxyaceic acid (PAA) for endoscopes washer. It is used for the automatized treatment of flexible endoscopes and other medical devices present in healt facilities or biologic labs.

H-Matic F

Concentrated high level disinfectant based on peracetic acid. For more info see page 20.

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DISINFECTANT based on peroxide and lactic acid

DESCRIPTION

LC-Act is a disinfectant solution. It contains 6% of Hydrogen Peroxide, alcoholic mix and Lactic Acid. Complete action spectrum. Effective on virus, bacteria, fungi, yeasts, spores.

COMPOSITION

6g Hydrogen Peroxide, 0,3 Lactic acid, Ethanol, Isopropanol, Citric Acid, Fosforic Acid, Demineralized water enough to do 100g.

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METHOD OF USE

Ready to use product. Product is used without dilution. Directly nebulize it on the medical devices.

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SPRAY DISINFECTANT hydroalcoholic solution

DESCRIPTION

Intermediate-level spray disinfectant, specific for disinfection and cleansing of medical devices. Hydroalcoholic formula. Available in wipes format.

COMPOSITION

Alcoholic mixture, o-phenylphenol, Didecyldimethylammonium chloride.

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LC-SAL METHOD OF USE

Nebulize with the appropriate dispenser directly on surfaces or supplies to be treated. Leave on for 1 minute for a bactericidal activity in clean condition. For further information consult the data sheet.

LC-SAL WIPES METHOD OF USE

Wipe directly on the surfaces and devices to disinfect using one or more wipes according to the area to be treated. Leave on at least for 5 minutes, dry if necessary. Do not rinse with water in order to get a residual biocidal effect on the surfaces. Do not use the same wipe more than once.

LC-Guan e LC-Guan Wipes

✓ Formula inodore
Scentless formula

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SPRAY DISINFECTANT solution

DESCRIPTION

Spray disinfectant solution with fungicidal and bactericidal activities, active on enveloped viruses. Specific for sanitization, decontamination and disinfection of non-invasive medical devices. Aqueous solution. Also available in wipes format.

COMPOSITION

Polyhexamethylene biguanide, Didcyldimethylammonium chloride.

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LC-SAL METHOD OF USE

Nebulize with appropriate spray dispenser directly on the surfaces or supplies to be treated. For further information consult the data sheet.

LC-SAL WIPES METHOD OF USE

Wipe directly on the surfaces and devices to disinfect using one or more wipes according to the area to be treated. Leave on at least for 5 minutes, dry if necessary. Do not rinse with water in order to get a residual biocidal effect on the surfaces. Do not use the same wipe more than once.

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Spray DISINFECTANT alcol based

DESCRIPTION

Septavir Alcol 75+ is an alcohol solution, ready to use, without rinse, with disinfection action, recommended for the disinfection of all the highly contaminated washable surfaces, like health facilities, hospitals, highly turnout business activities.

COMPOSITION

100g of product contains: Active ingredient Ethanol 75g, cocormulants, perfume: limonene, water enough to do 100g.

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SEPTAVIR ALCOL 75 METHOD OF USE

Nebulize 10ml of product for m² of clean surface, let it act for 2 minutes for an effectiveness on bacteria, for 5 minutes for a complete effectiveness (virus, fungi, yeasts and micobatteri). Let it dry. in case of surfaces that could be in contact with food, rinse with enough water.

Septavir Superfici

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- ✓ Nessun alone
No marks
- ✓ Senza solventi
Without solvents

- 📦 **Packaging**
 - Cod/SKU SSPMC/750
Fiacone da 750 ml/Bottle 750 ml
 - Cod/SKU SSPMC/1000
Fiacone da 1000 ml/Bottle 1000 ml
 - Cod/SKU SSPMC/5000
Tanica da 5000 ml/Jerrycan 5000 ml

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Spray DISINFECTANT based on Benzalkonium chloride

DESCRIPTION

SEPTAVIR SUPERFICI is a detergent solution, ready to use, based on Benzalkonium chloride. It removes germs and bacteria, significantly reducing the risk of their proliferation. Recommended for the daily disinfection of hard surfaces, walls, doors, at highly contamination risk too.

COMPOSITION

100 g of product contains: Active ingredient (benzalkonium chloride 0,45 g), surfactants, Co-formulated, perfume and water 100 g.

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SEPTAVIR SUPERFICI METHOD OF USE

Nebulize the product directly on the surface to be treated and let it act for 5 minutes for a bactericidal activity. For a complete activity on virus and fungi too let it act for 20 minutes, afterwards dry with clean cloth or paper. For surfaces that could be in contact with food rinse with water.

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DISINFECTANT based on Quaternary ammonium salts

DESCRIPTION

SEPTAVIR PAVIMENTI is a concentrated disinfectant solution based on Quaternary ammonium salts. Specific for cleaning and disinfection of floors and surfaces. The solution hasn't got corrosive characteristics thanks to the presence of antioxidants and it release pleasant scenting.

COMPOSITION

100 g of product contains: Active ingredient (Benzalkonium chloride 1,5 g) surfactants, phosphates, preservatives, perfume and water enough to do 100 g.

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SEPTAVIR PAVIMENTI METHOD OF USE

Hard surfaces disinfections: for a deep disinfection and sanitizing action in case of extreme dirtiness (EN 13697) dilute 25% (250 ml for every liter of wather). Let it act for 15 minutes and then rinse. For a disinfection action versus fungi and bacteria in accordance with the rules EN1276 - EN 1650: 7,5 (75ml for every liter of water). Let it act for 15 minutes and then rinse. For the daily hygiene of the floors: add 50ml of product in half bucket of water (4 liters). It doesn't need rinse.

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Surgical procedures may be carried out one after the other without the need for handwashing, provided that the handrubbing technique for surgical hand preparation is followed (images 1 to 17).

1
Put approximately 5ml (3 doses) of alcohol-based handrub in the palm of your left hand, using the elbow of your other arm to operate the dispenser

2
Dip the fingertips of your right hand in the handrub to decontaminate under the nails (5 seconds)

3
Images 3-7: Smear the handrub on the right forearm up to the elbow. Ensure that the whole skin area is covered by using circular movements around the forearm until the handrub has fully evaporated (10-15 seconds)

4
See legend for image 3

5
See legend for image 3

6
See legend for image 3

7
See legend for image 3

8
Put approximately 5ml (3 doses) of alcohol-based handrub in the palm of your right hand, using the elbow of your other arm to operate the dispenser

9
Dip the fingertips of your left hand in the handrub to decontaminate under the nails (5 seconds)

10
Smear the handrub on the left forearm up to the elbow. Ensure that the whole skin area is covered by using circular movements around the forearm until the handrub has fully evaporated (10-15 seconds)

11
Put approximately 5ml (3 doses) of alcohol-based handrub in the palm of your left hand, using the elbow of your other arm to operate the distributor. Rub both hands at the same time up to the wrists, and ensure that all the steps represented in Images 12-17 are followed (20-30 seconds)

12
Cover the whole surface of the hands up to the wrist with alcohol-based handrub, rubbing palm against palm with a rotating movement

13
Rub the back of the left hand, including the wrist, moving the right palm back and forth, and vice-versa

14
Rub palm against palm back and forth with fingers interlinked

15
Rub the back of the fingers by holding them in the palm of the other hand with a sideways back and forth movement

16
Rub the thumb of the left hand by rotating it in the clasped palm of the right hand and vice-versa

17
When the hands are dry, sterile surgical clothing and gloves can be donned

Riferimento: 2009 WHO Guidelines on Hand Hygiene in Health Care

The handrubbing technique for surgical hand preparation must be performed on perfectly clean, dry hands. On arrival in the operating theatre and after having donned theatre clothing (cap/hat/bonnet and mask), hands must be washed with soap and water. After the operation when removing gloves, hands must be rubbed with an alcohol-based formulation or washed with soap and water if any residual talc or biological fluids are present (e.g. the glove is punctured).

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DISINFECTANT alcohol-based

DESCRIPTION

Septavir Gel is a colourless alcoholic solution, in gel form, ready to use for disinfection of hand in absence of water. For every situation in which is necessary and suggested the sanification of hands for health workers and not only.

COMPOSITION

Ethyl alcohol, citrus fruit perfume.

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METHOD OF USE

Septavir gel is practical and ready to use. Doesn't need water and doesn't have to be rinsed after use. Antisepsis of hands. Rub 3 ml of products on dry hands for 30 seconds and let it dry. Repeat the operation a second time, for a total of 60 seconds.

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SANITIZING HAND SOAP based on lactic acid

DESCRIPTION

SEPTAVIR SOAP is a lactic acid based sanitizing detergent hand soap that physically removes the bacterial flora away the microbial flora of the skin. It's suitable for hand washing of medical staff, medical and dental clinics and all the facilities that need a high hands' hygiene rate.

COMPOSITION

Lactic acid, surfactants, lavender and olive oil essence.

METHOD OF USE

to cleanse, pour a small amount (5ml) of soap in the hollow of the hands and rub for at least 1 minute, rinse again. If necessary repeat the operation also 10 times/day.

Septavir Chlo Alcolico

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✓ Efficacia virucida completa
Complete efficacy virucides

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DISINFECTANT based on chloroxidine alcohol

DESCRIPTION

SEPTAVIR CHLO ALCOLICO is a ready to use solution, it contains chloroxidine digluconate 2% (w/w) and ethyl alcohol 70% (w/w), for the disinfection of the skin and hands surgical antiseptics. Glycerine contained in it helps to prevent the dryness of the skin, that could happen after the frequent usage of alcoholic solutions.

COMPOSITION

100 g of product contains: Active ingredients: chloroxidine digluconate 2g, ethyl alcohol 70g, co-formulants, water enough to do 100g.

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METHOD OF USE

For the disinfection of skin areas, put on the product directly on the skin or using a soaked cotton swab. Let it act for 1 minute for a virucidal activity (poliovirus, adenovirus, norovirus).

- Hands antiseptics: on dry hands, rub 3ml of product for 60 seconds and let it dry.
- Surgical hands antiseptics: put on the hands for 2 times 5 ml of product rubbing 2 minutes every time (total time 4 minutes).

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DISINFECTANT based on iodopovidone

DESCRIPTION

Antiseptic - antibacteria to use in the following cases:
 Disinfection of the skin. Antiseptic-microbiocide usage in
 medical clinics.

COMPOSITION

100 ml of solution contains: Polyvinylpyrrolidone iodine
 10g (equivalents to 1g of iodine) coformulated and
 demineralized water enough to 100 ml.

METHOD OF USE

Use it pure putting on directly on the interested area with
 a cotton baton or with a swab or gauze about 10ml for 1
 minute

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DISINFECTANT SOAP based on clorexidine

DESCRIPTION

It's a disinfectant liquid soap with bactericidal action
 in a ready to use soap solution based on clorexidine
 digluconate 4%.

COMPOSITION

100g of product contains: clorexidine digluconate 4g;
 perfume; co-formulated and demineralized water enough
 to do 100g.

METHOD OF USE

Use 5ml of product in contact with moistened hands for a
 minute. For the surgical preparation of the hands, extend
 the contact time to 2 minutes.

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REPAIR CREAM

DESCRIPTION

Based on: lenitive aloe vera; avocado oil that plays a sweetening and protective action versus the aggression of the weathering agents; karitè butter with emollient properties, useful in case of dermatitis, irritations and skin cracks; regenerating-acting urea.

COMPOSITION

AQUA, CETEARYL ALCOHOL, GLYCERIN, DICAPRYLYL ETHER, PERSEA GRATISSIMA OIL, UREA, ALOE BARBADENSIS LEAF JUICE, GLYCERYL STEARATE CITRATE, PROPYLENE GLYCOL, BUTYROSPERMUM PARKII BUTTER, CETEARETH-25, PARFUM, TOCOPHEROL, BENZYL ALCOHOL, SALICYL ACID, SORBIC ACID, DEHYDROACETIC ACID, PARFUM, XANTHAN GUM, DISODIUM EDTA, SODIUM HYDROXIDE, BENZYL SALYCILATE, D-LIMONENE, HEXYL CINNAMAL, LINALOOL, COUMARIN.

METHOD OF USE

massage evenly on the hands until completely absorbed.

MATERIALE ASSOCIATO
ASSOCIATED EQUIPMENT

**PROF**

- ✓ Formula concentrata
Concentrated formula
- ✓ Igienizza e deodora
Sanitizes and deodorizes

Packaging
 Cod/SKU PRFA/300
 Flacone da 300 ml/Bottle 300 ml
 Cod/SKU PRFC/300
 Flacone da 300 ml/Bottle 300 ml
 Cod/SKU PRFS/300
 Flacone da 300 ml/Bottle 300 ml
 Cod/SKU PRFV/300
 Flacone da 300 ml/Bottle 300 ml

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PROFESSIONAL SANITISING PERFUMER

DESCRIPTION

PROF® is a concentrated perfumer for professional usage and it's been formulated for deodorise and sanitise the highly turnout public spaces like waiting rooms, public transport, cinema, restaurants, ecc. Innovative formulation with Neutral pH, with Anti-static action, Able to removes bad smells, scent, clean, thoroughly sanitise, removing germs and bacteria and dust away every surfaces. The fragrances, accurately selected, persist over time.

COMPOSITION

(REG. CE 648/2004): (REG. CE 648/2004): Non-ionic surfactants % (w/w) < 5; Perfumes. Preservatives: DIDECYLDIMONIUM CHLORIDE, METHYLCHLOROISOTHIAZOLINONE / METHYLISOTHIAZOLINONE.

METHOD OF USE

Nebulize the product on a damp cloth for the cleaning of leather. Diluite 1ml of product in one liter of water for scent the floors, sanitise and dust away the surfaces. Before the usage, test the product on a little hidden area of the surfaces to treat.

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DETERGENT FOR GLUE

DESCRIPTION

LC-SOLV is a ready to use solvent for the removal of heavy dirty and glue residues.

COMPOSITION

100 g of product contains: solvent, excipients enough to do 100g.

METHOD OF USE

Put on with a sponge or a dry cloth, wait a few minutes and rinse. Put on directly or with a sponge, wait a few minutes and remove.

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PRODUCT INDEX in alphabetical order

■ Accessori 42	■ LC-Guan 28
■ Benzaclor 9	■ LC-Guan Wipes 28
■ Benzaclor P 9	■ LC-Sal 27
■ Crema mani riparatrice 38	■ LC-Sal Wipes 27
■ H-Matic ACT 23	■ LC-Solv 41
■ H-Matic AD 17	■ Neosterixidina 37
■ H-Matic Clean 18	■ Neozim 10
■ H-Matic DET 23	■ Neozim D 10
■ H-Matic F 20	■ Oxicer 11
■ H-Matic FI5 20	■ Pathoster® 12
■ H-Matic LUB 21	■ Pathoster® Matic 22
■ H-Matic NP 19	■ PROF 40
■ H-Matic NPC 19	■ Rugilux 15
■ H-Matic Rinse 21	■ Septavir Alcol 75+ 29
■ H-Matic SAD 17	■ Septavir Chlo Alcolico 36
■ H-Matic Zyme 18	■ Septavir Gel 34
■ H-Zim 0.5 7	■ Septavir Pavimenti 31
■ H-Zim P 8	■ Septavir SOAP 35
■ Iodopovidone 37	■ Septavir Superfici 30
■ LC-Act 25	■ Stericer 14

■ STRUMENTARIO

■ DISINFETTANTI

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■ MACCHINE

■ ANTISETTICI e DETERGENTI

A photograph of a modern, two-story industrial building with a white facade and blue accents. The building features large windows with blue frames and a prominent logo on the upper wall. The sky is blue with scattered white clouds.

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Cerichem Biopharm may change packaging and colors catalog items at any time and without prior notice. Any changes do not affect product effectiveness.

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Cerichem Biopharm S.r.l. - Viale Einaudi, snc (Z.I.) - 71042 Cerignola (FG)